



MG Scientific

8500 107th Street, Pleasant Prairie, WI 53158 • Phone: 262.947.7000 Fax: 262.947.7007



EPA UCMR5 Compliant

## Certified PFAS-Free Containers

Minimize contamination risks and maximize accuracy in your test results.

### EPA METHOD 1633

250 mL Natural HDPE Oblong LEAKPROOF, 45PP Cap,  
300 per Case

**B668-337**

### EPA METHOD 537.1

250 mL Natural HDPE Oblong LEAKPROOF, 45PP Cap,  
TRIS/TRIS HCL (1.25 g), 300 per Case

**B668-338**

### EPA METHOD 533

250 mL Natural HDPE Oblong LEAKPROOF, 45PP Cap,  
Ammonium Acetate (250 mg), 300 per Case

**B668-339**

[MGScientific.com/PFAS](http://MGScientific.com/PFAS)  
Robert McGregor  
262-237-9639